Become a SomaLogicTM Certified Site




Reliable Sample Insights

SomaScan[™] 7K Assay Kits detect the most proteins over a broad dynamic range, without sacrificing reproducibility, sensitivity, or specificity. Equipping your laboratory with kits provides:

- On-demand, on-site access to large-scale proteomics • Experimental control for high-throughput projects
- Easier workflow and less "hands on" lab time compared with other assays

Assay Specifics




Samples and controls per assay vary based on the kit type ²Typical results

Assay Workflow

omalogic

5





1FTE can run

340

samples/per week

With staggered runs, 2 FTE can run

680

samples/per week

3.5 hours

0.5 hours

Data Review and Upload

- Documentation (20 min)
- Data review and upload (10 min)

Hands On Steps Hands Off Steps Note that some hands-on and hands-off tasks run concurrently.



*Provided as an estimate. ADAT file-generation lead-time estimates are subject to fluctuate depending on demand.

O SomaLogic normalizes data, assures data

quality, and creates an ADAT file

O SomaLogic flags any data irregularities

ADAT File Generation

SL00000712 Rev 2: 2022-09

Become a SomaLogic Certified Site - Infographic

SomaLogic, SomaScan, SOMAmer, SomaSignal and associated logos are trademarks of SomaLogic Operating Co., Inc. and any third-party trademarks used herein are the property of their respective owners. © 2022 SomaLogic Operating Co., Inc. | 2945 Wilderness Pl, Boulder, CO 80301 | Ph 303 625 9000 | www.somalogic.com